




SIX MINUTES TO YOU COULD MEAN  
A LIFETIME TO HER.



The NEW Gammex 507 Rotational Therapy Phantom performs seven crucial quality control checks in just six minutes.

The level of care for your patients has reached new heights – and dimensions – with rotational therapy. However, it's difficult to ensure rotational accuracy with standard QA methodology, which is not designed to test the 360 degree interplay between all the systems. In just six minutes, the 507 Rotational Therapy Phantom performs quality checks on seven critical areas in rotational mode. As effective as it is efficient, the phantom and its software allows inspection of the MLC and gantry position synchrony, beam profile comparisons, constancy checks, and dose rate checks. When you're seeking a method to ensure accurate rotational therapy, and your patient's peace of mind...Think Gammex First.



For more information contact a  
GAMMEX representative today at  
**1 800 GAMMEX 1 (426 6391)**

**VISIT US AT  
AAPM BOOTH #1018**

